

## Questions & Answers

**Q.** How soon will I see results?

**A.** Some individuals will see immediate changes. But as with most exercise programs, the improvement becomes more noticeable after a period of time.

**Q.** Is the procedure safe?

**A.** Absolutely! The unit is battery powered and completely safe and easy to use. The only sensation that may be noticed is a slight pulsating feeling. As a precaution, those who have pacemakers will not want to use the Toning Unit below the chin area.

**Q.** How Long do the results last?

**A.** Like all exercise, the more you exercise your muscles on a regular basis, the better the results.

**Q.** How often should I use the system?

**A.** You can use it daily, but you will gain results by using it 2 to 3 times a week.

**Q.** What will happen to my face if I discontinue using the system?

**A.** Like all exercise programs, when you stop, the muscles will begin to lose the increased tone and over a period of time, the face will return to its former appearance.

**Q.** Will I know how to use the Toning Unit?

**A.** Included with each Toning Unit is a step-by-step guide along with 5 specialty products to enhance the toning procedure.

**Q.** Is there a warranty?

**A.** The Unit is warranted for 1 year against defects and workmanship. LYC Classic will repair any defective component for one full year for a small handling and shipping charge.

## LYC Classic Makers of the Facial Toning System And many other fine products



Visit us on the Internet at  
<http://www.lyconline.com>  
to learn about more of our fine,  
time-tested products.

**Distributed by:**



# LYC Classic

## Toning System

The Ultimate Home Facial Exerciser



**New! Now comes with a wonderful travel bag.**

**LYC CLASSIC**  
Stockton, CA, USA

## Can You Reverse The Aging Process?

The aging process has many visible indicators. The most obvious being our skin condition. No one likes to think of themselves as looking older. Beautiful, healthy-looking skin has been the desire of women through the ages. Today, more and more men are becoming concerned about their skin's health and appearance as well.


There are many factors that can contribute to your skin's appearance, such as poor dietary habits, every day pressures, heredity and environmental aggressions. Another major factor that contributes to the wrinkling and sagging of the skin is the gradual deterioration of muscle tone and a loss of resiliency in the connective tissues located between the facial muscles and the skin.

It is known that with age, muscles atrophy and weaken. As muscle fibers deteriorate, the muscle shrinks and the skin attached to that muscle begins to sag and wrinkle.

## How Does The LYC Classic Toning System Work?

The LYC Classic Toning Unit utilizes micro-current technology that sends a mild, pulsating current into the skin to help exercise and restore tone to facial muscles.

This process helps to diminish fine lines and increase circulation. It also firms, tightens and lifts your facial muscles to a more radiant and healthy appearance.



In addition, the micro-current helps to carry nutrients into the underlying layers of the skin. The result is a dramatically improved skin texture and a return to a glowing look of youth.

## Now You Can Have the Face Of The Future

For centuries Ancient China has produced many special forms of treatment that helps people's health and beauty. The LYC Classic Toning System utilizes a way of reaching the various acupuncture points in the face that have long been known in the Orient as "youth" points to retain a more healthy, youthful and glowing appearance.

Now in the privacy of your home and by spending 10-15 minutes a day, you can exercise and stimulate these points and produce beautiful results.

The LYC Classic Toner should be used to exercise the muscles around the chin, jaw, mouth, eyes, forehead, ear and throat. It can also be used in the back of the neck and shoulder area to help reduce the tension from stress that tightens the muscles.

